



Be Right™



Electrode Cleaning solution, RENOVO.N, for Clean Water Samples, 250 mL

Product #:

S16M001

Hazardous

GBP Price:
Available

Contact Us



Electrode cleaning solution for pH sensing element regular cleaning maintenance. Convenient for soft cleaning of the pH sensing element after use in clean water samples. 250 mL bottle.

Specifications

Description:	After measurements in aqueous solutions, cleaning with a mild solution such as RENOVO.N is sufficient.
Formula:	KATHON CG, DECONEX 11
Parameter:	Cleaning solution
Quantity:	250 mL
Shelf life:	18 months from production date
Solution Type:	Cleaning