



Be Right™



Hardness (Total & Calcium) Test Kit, Model HAC-DT

Product #: 2063900
GBP Price: Contact Us
Call for ship date

Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 10 - 4,000 mg/L. Contains approximately 100 of each (Total and Calcium) test.

Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration

Rugged titration device is ideal for field use, with no glass burets required

Specifications

| | |
|---------------------|--|
| Case style: | F |
| Measuring range: | 10 - 4,000 mg/L CaCO ₃ |
| Method: | Digital Titrator/EDTA |
| Number of tests: | 100 Ca +100 T |
| Parameter: | Hardness, Total & Calcium - as CaCO ₃ |
| Platform: | Digital Titrator |
| Shipping weight: | 4 |
| Smallest increment: | 0.1-10 |
| Test kit model: | HAC-DT |

What's in the box

Includes 6 reagents, digital titrator, delivery tubes, graduated cylinder, flask, manual and carrying case.