



Be Right™



Ammonium cuvette test 2.0-47.0 mg/L NH₄-N, 25 tests

Product #:
GBP Price:
Available

LCK303
Contact Us

Hazardous



LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Officially approved for consent limits. Environmentally friendly due to less chemicals.

Precise and reliable measurement values

Maximum safety for users, thanks to the closed system and low amounts of reagents

Convenient and error-free dosing of the reagents

Barcode label for automatic recognition in the photometer and Truecal calibration

Differentiation between tests and measuring ranges by means of colour-coding

Specifications

According to standard:	ISO 7150-1, DIN 38406 E5-1, UNI 11669:2017
EPA compliant:	N/A
Instrument:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000
Measuring range:	2 - 47 mg/L NH ₄ -N
Measuring range (2):	2.5 - 60 mg/L NH ₄
Method:	Indophenol Blue
Number of tests:	25
Parameter:	Ammonium
Platform:	LCK
Shelf life:	18 months from production date
Storage conditions:	2 - 8 °C (keep refrigerated)
Truecal:	Yes