



Be Right™



## COD cuvette test 15-150 mg/L O<sub>2</sub>, 25 tests

Product #:

LCK314

GBP Price:  
Available

Contact Us

Hazardous



### LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Officially approved for consent limits. Environmentally friendly due to less chemicals.

#### Precise and reliable measurement values

Maximum safety for users, thanks to the closed system and low amounts of reagents

Convenient and error-free dosing of the reagents

Barcode label for automatic recognition in the photometer and Truecal calibration

Differentiation between tests and measuring ranges by means of colour-coding

---

## Specifications

According to standard:	ISO 6060-1989, DIN 38409-H41-H44
Description:	Chemical Oxygen Demand
Digestion required:	Yes
EPA compliant:	N/A
Instrument:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000
Measuring range:	15 - 150 mg/L O <sub>2</sub>
Method:	Dichromate
Number of tests:	25
Parameter:	COD
Platform:	LCK
Shelf life:	23 months from production date
Storage conditions:	15 - 25 °C (protect from light)
Truecal:	Yes