



Be Right™



## Nitrate cuvette test 0.23-13.5 mg/L NO<sub>3</sub>-N, 25 tests

Product #:

LCK339

GBP Price:  
Available

Contact Us

Hazardous



### LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Officially approved for consent limits. Environmentally friendly due to less chemicals.

#### Precise and reliable measurement values

Maximum safety for users, thanks to the closed system and low amounts of reagents

Convenient and error-free dosing of the reagents

Barcode label for automatic recognition in the photometer and Truecal calibration

Differentiation between tests and measuring ranges by means of colour-coding

---

## Specifications

According to standard:	ISO 7890-1-2-1986, DIN 38405 D9-2
EPA compliant:	N/A
Instrument:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000
Measuring range:	0.23 - 13.5 mg/L NO <sub>3</sub> -N
Measuring range (2):	1.0 - 60.0 mg/L NO <sub>3</sub>
Method:	2.6-Dimethylphenol
Number of tests:	25
Parameter:	Nitrate
Platform:	LCK
Shelf life:	23 months from production date
Standard method:	EN 38405 D-2
Storage conditions:	15 - 25 °C
Truecal:	Yes