



Be Right™



LT20 BOD5 thermostat

Product #:

LTV073

GBP Price:

Contact Us

Ships within 1 week

Temperature control unit for 27 sample cells, adjustable to 20°C ± 1°C.

Specifications

Dimensions (H x W x D):	185 mm x 200 mm x 230 mm
Fuse:	T 315 mA L; 250 V (x 2)
Operating temperature range:	10 - 30 °C max. 80 % relative humidity, non-condensing
Power consumption:	65 VA
Power requirements (Hz):	50 Hz
Power requirements (Voltage):	230 V
Storage conditions:	-20 - 60 °C max. 80 % relative humidity, non-condensing
Weight:	4 kg