



Be Right™



Ammonium cuvette test, 2-47 mg/L, for AP3900 Laboratory Robot

Product #:

APC303

GBP Price:
Available

Contact Us

Hazardous



Hach APC cuvette tests: Minimum amounts of reagents, maximum safety.

100 cuvettes designed for use with the AP3900 Laboratory Robot.

Measuring range: 2 - 47 mg/L NH₄-N.

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in 4 packs with 25 cuvettes each. Officially approved for consent limits. Environmentally friendly due to less chemicals.

Precise and reliable measurement values

Maximum safety for users, thanks to the closed system, low amounts of reagents and included recycling/disposal service

Convenient and error-free dosing of the reagents

Barcode label for automatic test recognition in the photometer, including batch number, documented shelf life, Hach's Truecal functionality and COA via RFID technology

Differentiation between tests and measuring ranges by means of colour-coding

Specifications

According to standard:	ISO 7150-1, DIN 38406 E5-1
Measuring range:	2 - 47 mg/L NH ₄ -N
Measuring range (2):	2.5 - 60 mg/L NH ₄
Method:	Indophenol blue
Number of tests:	100
Parameter:	Ammonium
Platform:	APC
Shelf life:	17 months from production date
Storage conditions:	2 - 8 °C (keep refrigerated)